

logP Plugin

This manual gives you a walk-through on how to use the logP Plugin:

- Introduction
- Options
 - General Options
 - Method
 - Other options
 - Display Options
- References

Introduction

The logP Plugin calculates the octanol/water partition coefficient, which is used in QSAR analysis and rational drug design as a measure of molecular hydrophobicity. The calculation method is based on the publication of [Viswanadhan et al.](#) The logP of a molecule is composed of the increment of its atoms. However, the algorithm described in the paper was modified at several points:

1. Many atomic types were redefined to accommodate electron delocalization and contributions of ionic forms were added.
2. The logP of zwitterions are calculated from their logD value at their isoelectric point.
3. The effect of hydrogen bonds on the logP is considered if there is a chance to form a six membered ring between suitable donor and acceptor atoms.
4. New atom types were introduced especially for sulfur, carbon, nitrogen and metal atoms.

To find details on logP calculation, see the [following page](#).

The result of the calculation appears in a new window, either in a MarvinView (for a [2D view](#)) window or in a MarvinSpace (for a [3D view](#)) window.




Fig. 1 Log P result window with atomic increments displayed in MarvinView




Fig. 2 Log P result window with atomic increments displayed in MarvinSpace

Options

General Options

Different general options can be set in the *logP Options* window:

Method

This option defines the method for calculating log P . These can be:

- **Consensus:** this method uses a consensus model built on the [ChemAxon](#) and [Klopman et al.](#) models and the [PhysProp](#) database.
- **ChemAxon:** this method is based on ChemAxon's own log P model, which is based on the VG method (derived from [Viswanadhan et al.](#)). Read more about it [here](#).

- **User defined:** if a training set of structures and corresponding experimental $\log P$ values is available, it can be used as a database for $\log P$ calculations. See [the manual page](#) on creating such training sets.

i The relation of ChemAxon's logP method to the AlogP and ClogP methods

The Consensus logP method is a unique, in-house developed logP model based on the methods listed above. Comparison-wise, our Consensus logP method is similar (but not identical) to the ClogP method, while our ChemAxon logP method is similar (but not identical) to the AlogP method.

Other options

These options are for refining the $\log P$ calculation.

- **Training ID:** if the **User defined** method is selected, this dropdown list becomes active. All created training sets are listed here. Choose the one to apply for the calculation.
- **Electrolyte concentration**
 - **Cl⁻ concentration:** can be set between 0.1 and 0.25 mol/L.
 - **Na⁺ K⁺ concentration:** can be set between 0.1 and 0.25 mol/L.
- **Consider tautomerization/resonance:** the $\log P$ of the major tautomer will be calculated if enabled.




Fig. 3 $\log P$ Options window with the General Options panel

Display Options

Different display options can be set in the *logP Options* window:

- **Decimal places:** setting the number of decimal places for the precision of the result value.
- **Show value:** sets which $\log P$ values be displayed in the result window. It can be:
 - **Increments:** calculates and shows the increments for atoms one-by-one.
 - **$\log P$:** shows just the value of the $\log P$.
- **Display in MarvinSpace:** the result window opens as [3D MarvinSpace viewer](#). If unchecked, the results will be shown in a [2D MarvinView panel](#).




Fig. 4 log P Options window with the Display Options panel

References

1. Viswanadhan, V. N.; Ghose, A. K.; Revankar, G. R.; Robins, R. K., *J. Chem. Inf. Comput. Sci.*, **1989**, *29*, 163-172; doi
2. Klopman, G.; Li, Ju-Yun.; Wang, S.; Dimayuga, M.: *J.Chem.Inf.Comput.Sci.*, **1994**, *34*, 752; doi
3. PHYSPROP[©] database
4. Csizmadia, F; Tsantili-Kakoulidou, A.; Pander, I.; Darvas, F., *J. Pharm. Sci.*, **1997**, *86*, 865-871; doi